


Product Datasheet

Mouse IL25 protein (orb80103)

Description	IL25 Mouse Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain containin...
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	mSF20
Preservatives	Mouse SF20 was lyophilized from 0.2 um filtered concentrated (1.0mg/ml) solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MVSEPTTVPF DVRPGGVVHS FSQDVGPNGK FTCTFTYASQ GGTNEQWQMS LGTSEDSQHF TCTIWRPQ GK SYLYFTQFKA ELRGAEIEYA MAYSKAAFER ESDVPLKSEE FEVTKTAVSH RPGAFKAELS KLVIVAKAAR SEL.
Purity	> 97.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized Mouse SF20 in sterile water not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse IL25 protein