


Product Datasheet

Human Trefoil Factor-3 protein (orb80106)

Description	TFF3 Human Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain containin...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	TFF3
Preservatives	The TFF3 protein was lyophilized from 0.4 um filtered solution at concentration of 0.5mg/mL containing 20 mM Tris pH-7.5 and 150mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MKHHHHHHAS EEYVGLSANQ CAVPAKDRVD CGYPHVTPE CNNRGCCFDS RIPGVPWCFK PLQEAECTF.
Purity	> 95.0% as by SDS-PAGE.
Solubility (25°C)	Add deionized water to working concentration approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by appropriate sterile filter before using it in the cell culture.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Trefoil Factor-3 protein.