


Product Datasheet

Human TNF alpha, Mutant protein (orb80111)

Description	Human TNF alpha, Mutant protein
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Human TNF alpha, Mutant
Preservatives	The protein was lyophilized after extensive dialysis against 0.5x PBS pH -7.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MRKRKPVAHV VANPQAEQQL QWLNRRANAL LANGVELRDNQLVVPSEGLY LIYSQVLFGK QGCPSTHVLL THTISRIAVS YQTKVNLLSA IKSPCQRETP EGAEAKPWYE PIYLGGVFQL EKGDRLSAEI NRPDYLDFAE SGQVYFGIIAF.
Purity	> 95.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized Tumor Necrosis Factor-alpha Variant in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human TNF alpha, Mu...