

Product Datasheet

Human B-Cell Maturation Antigen protein (orb80112)

Description	TNFRSF17 Human Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain conta...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	TNFRSF17
Preservatives	1mg of TNFRSF17 Human contain 20 mM sodium phosphate buffer, pH-7.4, and 130mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	AGQCSQNEYF DSSLHACIPC QLRCSSNTPP LTCQRYCNAS VTNSVKGTNA.
Purity	> 98.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized TNFRSF17 in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.