

## Product Datasheet

### Rat MIP2 protein (orb80133)

<b>Description</b>	GRO-Beta Rat Recombinant also caled Rat MIP-2 produced in E.coli is a single, non-glycosylated, poly...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FA, HPLC, SDS-PAGE
<b>Target</b>	rGRO b
<b>Preservatives</b>	The protein was lyophilized with no additives.
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	> 95.0% as determined by SDS-PAGE.
<b>Solubility (25°C)</b>	It is recommended to reconstitute the lyophilized GRO-beta Rat in sterile H <sub>2</sub> O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
<b>Expiration Date</b>	6 months from date of receipt.