

Product Datasheet

Viral MIP2 protein (orb80134)

Description	GRO-Beta Viral Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain containing...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	vGRO b
Preservatives	Lyophilized from 1mg/ml solution containing 20 mM Phosphate Buffer, pH 7.4 0.15M NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	LGASWHRPDK CCLGYQKRPL PQVLLSSWYP TSQLCSKPGV IFLTKRGRQVCADKSKDWVK KLMQQLPVTA.
Purity	> 97.0% as determined by RP-HPLC and SDS-PAGE analysis
Solubility (25°C)	It is recommended to reconstitute the lyophilized CXCL2 Viral in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.