

Product Datasheet

Human Lymphotactin protein (orb80139)

Description	Lymphotactin Human Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain co...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Lymphotactin
Preservatives	The XCL1 was lyophilized from concentrated (1mg/ml) solution in water containing no additives.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	The sequence of the first five N-terminal amino acids was determined and was found to be Gln-Ser-Glu-Val-Ser.
Purity	> 99.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized Lymphotactin in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.