

Product Datasheet

ND2 antibody (orb801416)

Description	ND2 antibody
Species/Host	Rabbit
Reactivity	Worm
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Ascaris suum (Pig roundworm) (Ascaris lumbricoides) ND2 protein
Target	ND2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24877
Expiration Date	12 months from date of receipt.