

Product Datasheet

Human TARC protein (orb80164)

Description	Human TARC protein
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Human TARC
Preservatives	The protein was lyophilized from concentrated (0.5mg/ml) solution containing 20 mM PBS 150mM NaCl pH-7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	ARGTNVGRECCLEYFKGAIPRLKLTWYQTSSEDCSRDAIVFVTVQGRAICSDPNN K RVKNAV KYLQSLERS.
Purity	> 97.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized CCL17 in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.