


Product Datasheet

Human TECK protein (orb80165)

Description	TECK Human Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain containin...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	TECK
Preservatives	Filtered (0.2 um) and lyophilized from concentrated (1mg/ml) solution in 20 mM PB, pH 7.4, 150mM NaCl.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	QGVFEDCCLA YHYPIGWAVL RRAWTYRIQE VSGSCNLPAA IFYLPKRHRK VCGNPKSREV QRAMKLLDAR NKVFAKLHHN MQTFQAGPHA VKKLSSGNSK LSSSKFSNPI SSSRKNVSL ISANSGL.
Purity	> 97.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized TECK in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human TECK
protein