

Product Datasheet

F14F9.5 antibody (orb802184)

Description	F14F9.5 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Caenorhabditis elegans F14F9.5 protein
Target	F14F9.5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9GUC9
Expiration Date	12 months from date of receipt.