

Product Datasheet

TFIP11 antibody (orb802650)

Description	TFIP11 antibody
Species/Host	Rabbit
Reactivity	Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey) TFIP11 protein
Target	TFIP11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A1XD95
Expiration Date	12 months from date of receipt.