

Product Datasheet

Y75B8A.31 antibody (orb802833)

Description	Y75B8A.31 antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Caenorhabditis elegans Y75B8A.31 protein
Target	Y75B8A.31
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9XW78
Expiration Date	12 months from date of receipt.