

## Product Datasheet

### K02A2.5 antibody (orb803158)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | K02A2.5 antibody                                                 |
| <b>Species/Host</b>        | Rabbit                                                           |
| <b>Reactivity</b>          | Other                                                            |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                        |
| <b>Immunogen</b>           | Recombinant Caenorhabditis elegans<br>K02A2.5 protein            |
| <b>Target</b>              | K02A2.5                                                          |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid<br>repeated freeze. |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG                                                              |
| <b>Clonality</b>           | Polyclonal                                                       |
| <b>Uniprot ID</b>          | <b>Q09401</b>                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |