

Product Datasheet

5-hydroxytryptamine receptor 1 antibody (orb803962)

Description	5-hydroxytryptamine receptor 1 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) 5-HT7 protein
Target	5-HT7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20905
Expiration Date	12 months from date of receipt.