

Product Datasheet

GF24309 antibody (orb804574)

Description	GF24309 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila ananassae (Fruit fly) GF24309 protein
Target	GF24309
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B3M6C8
Expiration Date	12 months from date of receipt.