

Product Datasheet

E(spl)m3-HLH antibody (orb805094)

Description	E(spl)m3-HLH antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) E(spl)m3-HLH protein
Target	E(spl)m3-HLH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01068
Expiration Date	12 months from date of receipt.