

Product Datasheet

Hermansky-Pudlak syndrome 5 protein homolog antibody (orb805646)

Description	Hermansky-Pudlak syndrome 5 protein homolog antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) p protein
Target	p
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9VHN9
Expiration Date	12 months from date of receipt.