

Product Datasheet

Hsc70-4 antibody (orb805910)

Description	Hsc70-4 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) Hsc70-4 protein
Target	Hsc70-4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11147
Expiration Date	12 months from date of receipt.