

Product Datasheet

IntS1 antibody (orb806012)

Description	IntS1 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) IntS1 protein
Target	IntS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9W1C5
Expiration Date	12 months from date of receipt.