

Product Datasheet

SCR antibody (orb807594)

Description	SCR antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster (Fruit fly) Scr protein
Target	Scr
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09077
Expiration Date	12 months from date of receipt.