

Product Datasheet

ARBH antibody (orb809748)

Description	ARBH antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Dickeya chrysanthemi (Pectobacterium chrysanthemi) (Erwinia chrysanthemi) arbH protein
Target	arbH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26217
Expiration Date	12 months from date of receipt.