

Product Datasheet

ARCA antibody (orb809750)

Description	ARCA antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) arcA protein
Target	arcA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7LMS3
Expiration Date	12 months from date of receipt.