

Product Datasheet

CBPM antibody (orb811063)

Description	CBPM antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks) cbpM protein
Target	cbpM
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B1IV98
Expiration Date	12 months from date of receipt.