

Product Datasheet

CSPA antibody (orb812263)

Description	CSPA antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti) cspA protein
Target	cspA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Z3S6
Expiration Date	12 months from date of receipt.