

Product Datasheet

CYNS antibody (orb812649)

Description	CYNS antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks) cynS protein
Target	cynS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B1J0T9
Expiration Date	12 months from date of receipt.