

Product Datasheet

DNA polymerase antibody (orb813754)

Description	DNA polymerase antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Bacillus phage M2 (Bacteriophage M2) G protein
Target	G
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19894
Expiration Date	12 months from date of receipt.