

Product Datasheet

tuf antibody (orb814254)

Description	tuf antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mycoplasma hominis (strain ATCC 23114 / NBRC 14850 / NCTC 10111 / PG21) tuf protein
Target	tuf
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22679
Expiration Date	12 months from date of receipt.