

Product Datasheet

HCAF antibody (orb817444)

Description	HCAF antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O139:H28 (strain E24377A / ETEC) hcaF protein
Target	hcaF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A7ZPY2
Expiration Date	12 months from date of receipt.