

Product Datasheet

INPP5A antibody (orb819069)

Description	INPP5A antibody
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Canis lupus familiaris (Dog) (Canis familiaris) INPP5A protein
Target	INPP5A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q29467
Expiration Date	12 months from date of receipt.