

Product Datasheet

NHA2 antibody (orb822853)

Description	NHA2 antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rhodococcus rhodochrous Low-molecular weight cobalt-containing nitrile hydratase subunit alpha protein
Target	Low-molecular weight cobalt-containing nitrile hydratase subunit alpha
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P29378
Expiration Date	12 months from date of receipt.