

Product Datasheet

P2OX antibody (orb823660)

Description	P2OX antibody
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Trametes pubescens (White-rot fungus) p2ox protein
Target	p2ox
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5G234
Expiration Date	12 months from date of receipt.