

Product Datasheet

PHHY antibody (orb823789)

Description	PHHY antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Trichosporon cutaneum (Yeast) (Oidium cutaneum) Phenol hydroxylase protein
Target	Phenol hydroxylase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15245
Expiration Date	12 months from date of receipt.