


Product Datasheet

Bacterial amiB2 protein (orb8243)

Description	Recombinant of bacterial amiB2 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag info: N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged Expression Region: 1-462aa Sequence Info: Full Length Glycerol content: 0.5 Relevance: A monocarboxylic acid amide + H ₂ O = a monocarboxylate + NH ₃ . Protein Description: Full Length
Protein Sequence	MDPTDLAFAGAAAQARMLADGALTAPMLLEVYLQRIERLDShLRAYRVVQFDRA RAEAEAAQRLDAGERLPLLGVPIAIKDDVDIAGEVTTYGSAGHGPAATSDAEVV RRLRAAGAVIIGKTNVPELMIMPFTESLAFGATRNPWCLNRTPGGSSGGSSAAAVA AGLAPVALGSDGGGSIRIPCTWCGLFGLKPQRDRISLEPHDGAWQGLSVNGPIA RSVMDAALLLDATTTVPGPEGEFVAAAAARQPGRRLIALSTRVPTPLPVRGKQEL AAVHQAGALLRDLGHDVVVRDPDYPASTYANYLPRFFRGISDDADAQAHPDRL EARTRAIARLGSSFFSDRRMAALRAAEVVLSSRIQSIFDDVDVVVTPGAATGPSRIG AYQRRGAVSTLLLVVQRPVYFQVWNLTGQPAAVVPWDFDGDGLPMSVQLVGR PYDEATLLALAAQIESARPWARRPSVS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	69.1 kDa
Uniprot ID	P9WQ97
Expiration Date	6 months from date of receipt.



Discontinuous
SDS-PAGE
(reduced) with
5%...