

Product Datasheet

PRVT antibody (orb824322)

Description	PRVT antibody
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Gallus gallus (Chicken) Parvalbumin, thymic protein
Target	Parvalbumin, thymic
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19753
Expiration Date	12 months from date of receipt.