

Product Datasheet

PP antibody (orb824403)

Description	PP antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia phage lambda (Bacteriophage lambda) Serine/threonine-protein phosphatase protein
Target	Serine/threonine-protein phosphatase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03772
Expiration Date	12 months from date of receipt.