


Product Datasheet

Human GCSF protein (orb82613)

Description	GCSF is a pleiotropic cytokine best known for its specific effects on the proliferation, differentia...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	GCSF
Preservatives	Lyophilized from a 0.2µm filtered concentrated (1mg/ml) solution in 10mM sodium acetate buffer pH 4.0. Lyophilized from a 0.2µm filtered concentrated (1mg/ml) solution in 10mM sodium acetate buffer pH 4.0.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	TPLGPASSLP QSFLKCLEQ VRKIQGDGAA LQEKLCATYK LCHPEELVLL GHSLGIPWAP LSSCPSQALQ LAGCLSQLHS GLFLYQGLLQ ALEGISPELG PTLDTLQLDV ADFATTIWQQ MEELGMAPAL QPTQGAMPAF ASAFQRRAGG VLVASHLQSF LEVSYRVLRLH LAQP
Purity	> 98% pure by SDS-PAGE and HPLC analyses.
MW	18.7 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human GCSF
protein