

Product Datasheet

rplQ antibody (orb826285)

Description	rplQ antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O7:K1 (strain IAI39 / ExPEC) rplQ protein
Target	rplQ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7NLL4
Expiration Date	12 months from date of receipt.