

Product Datasheet

rpsB antibody (orb827607)

Description	rpsB antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O8 (strain IA11) rpsB protein
Target	rpsB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7M1A9
Expiration Date	12 months from date of receipt.