

Product Datasheet

rpsT antibody (orb828550)

Description	rpsT antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC) rpsT protein
Target	rpsT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7N7P7
Expiration Date	12 months from date of receipt.