

Product Datasheet

YCHJ antibody (orb829271)

Description	YCHJ antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O139:H28 (strain E24377A / ETEC) ychJ protein
Target	ychJ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A7ZL06
Expiration Date	12 months from date of receipt.