

Product Datasheet

YMGJ antibody (orb830466)

Description	YMGJ antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O6:K15:H31 (strain 536 / UPEC) ymgJ protein
Target	ymgJ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q0TIK5
Expiration Date	12 months from date of receipt.