

Product Datasheet

YCR4 antibody (orb830669)

Description	YCR4 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia vulneris Uncharacterized 11.2 kDa protein in crtE 3'region protein
Target	Uncharacterized 11.2 kDa protein in crtE 3'region
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01333
Expiration Date	12 months from date of receipt.