

Product Datasheet

YCDZ antibody (orb830948)

Description	YCDZ antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) ycdZ protein
Target	ycdZ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O54290
Expiration Date	12 months from date of receipt.