

Product Datasheet

YFBU antibody (orb830983)

Description	YFBU antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O45:K1 (strain S88 / ExPEC) yfbU protein
Target	yfbU
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7MG58
Expiration Date	12 months from date of receipt.