

Product Datasheet

SMG antibody (orb832240)

Description	SMG antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) smg protein
Target	smg
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8FD15
Expiration Date	12 months from date of receipt.