

Product Datasheet

SOD antibody (orb832754)

Description	SOD antibody
Species/Host	Rabbit
Reactivity	Parasite
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Schistosoma mansoni (Blood fluke) SOD protein
Target	SOD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01137
Expiration Date	12 months from date of receipt.