

Product Datasheet

sodA antibody (orb832756)

Description	sodA antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli (strain K12) sodA protein (2-206aa)
Target	sodA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00448
Expiration Date	12 months from date of receipt.