

Product Datasheet

TATE antibody (orb832944)

Description	TATE antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280) tatE protein
Target	tatE
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	D2TMS1
Expiration Date	12 months from date of receipt.