

Product Datasheet

TATD antibody (orb833277)

Description	TATD antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Citrobacter rodentium (strain ICC168) (Citrobacter freundii biotype 4280) tatD protein
Target	tatD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	D2TUZ4
Expiration Date	12 months from date of receipt.